

## CROSS REACTIVITY - PCP

Various potential interfering substances in a human urine matrix were tested for cross-reactivity with the SYNCHRON Systems PCP assay. The following table summarizes the results obtained at the concentrations tested for each potential cross-reactant.<sup>1</sup>

**Table 3 Cross Reactivity<sup>2</sup>**

COMPOUND	CONCENTRATION (µg/mL)	EFFECT
Phencyclidine (cutoff)	0.025	Positive
Acetaminophen	1000	Negative
Acetylsalicylic Acid	1000	Negative
Amobarbital	1000	Negative
d-amphetamine	1000	Negative
Benzoylcegonine	1000	Negative
Brompheniramine	50	Negative
Chlorophiramine	50	Negative
Chlorpromazine	100	Negative
Dextromethorphan	1000	Negative
Diphenhydramine	100	Negative
EMDP (Methadone Metabolite) <sup>3</sup>	100	Negative
Imipramine	500	Negative
Ketamine	100	Negative
Meperidine	50	Negative
Methadone	1000	Negative
Methaqualone	100	Negative
Morphine	200	Negative
Naltrexone	10	Negative
Norpropoxyphene	100	Negative
Orphenadrine	200	Negative

---

COMPOUND	CONCENTRATION (µg/mL)	EFFECT
Oxazepam	1000	Negative
Phenobarbital	1000	Negative
Promethazine	100	Negative
Propoxyphene	1000	Negative
Thioridazine	80	Negative
Triprolidine	10	Negative

<sup>1</sup> It is possible that other substances and/or factors (e.g. technical or procedural) not listed above may interfere with the test and cause false results.

<sup>2</sup> Data shown was collected using SYNCHRON CX Systems. Equivalency between SYNCHRON LX Systems has been established by Deming regression analysis to SYNCHRON CX Systems.

<sup>3</sup> EMDP: 2-ethyl-5-methyl-3,3-diphenylpyrroline.